



# Pennsylvania Mobile

VETERINARY ULTRASOUND SERVICES



## HELPFUL REMINDERS FROM PENNSYLVANIA MOBILE VETERINARY ULTRASOUND SERVICES

*Our goal is to provide an efficient, seamless service with minimal disruption to the everyday flow of your hospital. Please review the below procedures for scheduling, preparing, and assisting with Pennsylvania Mobile VUS appointments to ensure a smooth service.*

### BEFORE YOUR APPOINTMENT:

- ✔ Call PA Mobile for appointment scheduling at **(443) 794-4700**
- ✔ Please complete the current submission form. This is included in the most recent PA Mobile marketing email. Please contact us and we will email the form if you do not have the current version.
- ✔ Patients should be dropped off by 9am the day of the procedure, having had **NO FOOD** after a normal meal 6pm the evening prior. Water is permitted.
- ✔ The morning of the procedure, your clinic should call PA Mobile to obtain a 2 hour window for the appointment. The sonographer will also call 20-30 minutes before arrival. This phone call is for the front staff to communicate to the doctor/technician overseeing the patient, to ensure that the patient is **shaved, sedated, and all case history has been submitted**. Please feel free to contact us with any sedation recommendations.

### DURING YOUR APPOINTMENT:

- ✔ If the patient is prepped and ready to be scanned, with the assistance from 1 staff member, the scan should take approximately 15 minutes. Additional cavities and sampling take additional time.

### AFTER YOUR APPOINTMENT:

- ✔ Reports are available within 24 hours on the Sonopath Reports site ([reports.sonopath.com](https://reports.sonopath.com)). An auto generated email will be sent when it is available for viewing/downloading. If it has been over 24 hours and your report is not accessible on the site, please contact PA Mobile immediately.

Pennsylvania Mobile Veterinary Ultrasound Services, LLC  
443.794.4700  
[rkjakum@gmail.com](mailto:rkjakum@gmail.com)  
[www.pamobilevet.com](http://www.pamobilevet.com)